

Simple Lab Clave
KYR-23

Horizontal laboratory autoclave having durability with simple operation and personal use.



KYR-23

Horizontal and small capacity table top autoclave

Compact and space saving model

Suitable for a wide range of applications

Chemical plants, Bio-Labs, Water & Air pollution Lab, Fisheries, Food QC Mass kitchens, etc.

“Mecha minimum ism” fits a minimum budget

No performance compromises even on the minimum budget.

MODEL	KYR-23
Type of placement	Table-top
Chamber size/volume	φ230×H460mm / 22L
Outer dimension/weight	W455×D520×H370mm / 29kg
Temperature(Max. press.)	110°C~125°C (0.160MPa)
Sterilization time	1 minite 00 seconds~99 minites 59seconds
Safety device	Water lack prevention · Over pressure protection · Electric leakage breaker · Pressure safety valve
Power demand	AC220/230/240V, 1 phase, 8A
Accessories	Stainless tray L(large) 1 pc. / Stainless tray S(small) 1pc.